


1. Remove the Myriad Clinical Pack Overwrap bag and the mailing label (found beneath the cardboard tray).

2. Place completed **TEST REQUEST FORM** with **PATHOLOGY REPORT** in the test kit. Please **do not** write the patient's name on the test kit.



3. Close the box and place it into the Myriad Clinical Pack Overwrap bag. Affix the enclosed pre-paid mailing label onto the front of the bag.*

4. If you don't have a regularly scheduled FedEx pickup, call 800-GO-FEDEX (800-463-3339) to request one. Press "0" and specify that you are shipping with a prepaid air bill.



Once we receive your test, Myriad will automatically replenish your kit supply. If you would like to order additional tests or have any questions about your specimen submission, please call Myriad myChoice CDx Customer Service at 1-877-283-6709.

*It is the shipper's responsibility to ensure that the package containing a diagnostic specimen conforms to FedEx guidelines along with all applicable local, state, federal and international regulations.

